

COLOUR CODE CHART

Restorative (Highspeed handpieces)

	Room #1	Room #2	Room #3	Room #4
Bottom Band (Room Designation)	No band	Black - ●	Gray - ●	Brown - ●
	1 - ●	1 - ●	1 - ●	1 - ●
Top Band (Use Order)	2 - ●	2 - ●	2 - ●	2 - ●
	3 - ●	3 - ●	3 - ●	3 - ●
	4 - ●	4 - ●	4 - ●	4 - ●
	4 - ●	4 - ●	4 - ●	4 - ●

Hygiene (Slowspeed attachments)

	Room #5	Room #6	Room #7	Room #8
Bottom Band (Room Designation)	Pink - ●	Yellow - ●	Mauve - ●	White- ●
	1 - ●	1 - ●	1 - ●	1 - ●
Top Band (Use Order)	2 - ●	2 - ●	2 - ●	2 - ●
	3 - ●	3 - ●	3 - ●	3 - ●
	4 - ●	4 - ●	4 - ●	4 - ●
	4 - ●	4 - ●	4 - ●	4 - ●

Sticker on the outside of each unit will allow staff to easily find correct room.



Colour Coding & Storage System

Optimize service life and reduce maintenance costs by ensuring even wear and use of all the handpieces in your practice.



Instructions and Colour Code Guidelines

Installation to wall:

- 1) Place The Handpiece Closet in the exact position where you wish to have it secured to wall or cupboard.
- 2) Place a level on top of The Handpiece Closet. With a pencil, mark the pre drilled holes to the wall or cupboard in order to locate exact location for drilling.
- 3) Using a ¼" drill bit, drill the locations of each mark that you just made.
- 4) If installing to drywall, use a small hammer to gently tap anchor into each hole.
- 5) Place The Handpiece Closet in exact position and apply screw to each hole for secure installation.

How The Handpiece Closet works:

- 1) Place coloured ID bands on each handpiece.
- 2) TOP coloured ID band for each handpiece (Blue, Orange, Green, Red) are to match the coloured sticker on each tray.
- 3) BOTTOM coloured ID band for each handpiece (NONE, Black, Gray, Brown, Pink, Yellow, Mauve, White) are to match the room that the handpiece is used in, and the round sticker on the outside of each Handpiece Closet.



- 4) Restorative - If using slowspeed attachment with highspeed, have exact colour ID bands for each so the pair always stay together.



- 5) Hygiene - Place ID bands on neck of head only. No need to have ID bands on contra angle as well. When sterilizing, separate head and drive shaft from contra angle, bag attachments and then sterilize.
- 6) To get started, place handpiece with proper coloured ID bands onto matching coloured tray in matching coloured room.
- 7) Starting from #1 Blue, use each handpiece in sequence. Never use handpieces out of numerical or coloured coded sequence.



Example: If you finish your day and all is left in The Handpiece Closet is #4 Red then that is the handpiece that you start off with the very next day.

If The Handpiece Closet is empty at any time, that indicates to start again from the top with #1 Blue. Refill The Handpiece Closet and start again, working your way through the proper colours, 1 through 4.

The purpose of this system is to have ALL handpieces used and sterilized the exact same amount of times as every other handpiece. NEVER USE COLOURS AND NUMBERS OUT OF SEQUENCE.

This system is guaranteed to reduce the cost of handpiece repair, saving your office thousands of dollars during years of practice.

www.thehandpiececloset.ca Toll free: 1-866-591-7770